

Non-Invasive Glucose Device



Group 18

Vithooban Thavapalan, Elia Tohme, Nicolas Freundler, Simon Romanski



What is diabetes ?

Diabetes mellitus :

Bad glucose uptake by cells

No insulin (type I)

No sensitivity to insulin (type II)

Consequence: Hyperglycaemia

Impact:

442.000.000 people in 2014

7th cause of death in 2030




Existing technology

Stripes



Continuous Glucose Monitoring (CGM)



→ All existing approaches are **invasive**



Our Solution


A non-invasive, wearable & discrete glucose monitoring device



Detection Principle: **NIR**

Absorption & Scattering

Wavelength & molecule dependent




$$I = I_0 \exp(-\mu_{eff} L)$$



Detection Principle: **NIR**



Glucose absorption spectrum



Used wavelength
940 nm



Detection Principle: **NIR**



http://ieeexplore.ieee.org/xpls/icp.jsp?arnumber=6777023&tag=1#ref_4




Our Device






Our Device






Concept of our device







Concept of our device






Concept of our device






Concept of our device






Concept of our device






Circuit Concept




Measurement






Circuit Concept




Measurement

Amplifier + Filter





Circuit Concept





Circuit Concept

D1: LED - 940nm (LI1IZ940)


D2: PV-Diode (PIN-3CDP)

R1 = 1Ω

R2 = $500M\Omega$

C1 = 16pf

C2 = 100nF






Our Device

Pulse	500	ms
Pulse every	15	min
Usage per day	48	s
LED forward current	1000	mA
Bluetooth Peak Current	20.6	mA
Op-Amp	<1	mA
Total	< 1025	mA
Battery Capacity	48	mAh
Recharge after	60	days



Our Device




Benefits

Non-Invasive



Our Device



Benefits

Non-Invasive

Convenient



Our Device



Benefits

Non-Invasive

Convenient

Discrete





Manufacturing - Flexibility

LED




Thickness : 700 μm



Interconnections

BENDING



**WITH RIGID ISLAND
UNDERNEATH**

ISLANDS : Relieve stress onto existing rigid components




Manufacturing - Stretchability

Bridge-Island Technique:

1. Adding a rigid island under LED
2. Keeping serpentine bridge

Expands the stretchability even further





Manufacturing - Materials



Substrate:

PDMS

Rigid Island:

Raisin Epoxy

Thin metal layer for interconnection:

Platinum or gold


Encapsulation:

PDMS



Manufacturing – Process flow


Manufacturing of
2 parts






Manufacturing – Process flow

First Part




1. Create the matrix



2. Cast PDMS on matrix



3. Bake PDMS
4h | 75° C



4. Flip the PDMS over




5. Pick and Place rigid
components




Manufacturing – Process flow


Second Part



1. Place islands




2. Cast PDMS on islands



3. Flip the PDMS over



4. Cast Gold/Platinum on layers



5. Shadow Mask Gold



Thank you